

ES cell genomic DNA purification in 96 well format

The protocol below is a quick purification of ES cell genomic DNA. The protocol uses a 3hr digestion of the ES cells, followed by isopropanol precipitation of the DNA.

Reagents (All reagents are from Sigma)

GenLysis Buffer

Stock concentration	volume (ml)	Final concentration
1M TrizmaHCl pH8	1ml	20mM Tris/HCL pH8
0.1M EDTA pH8	0.5ml	1mM EDTA
5M NaCl	0.05ml	5mMNacl
20% SDS	2.5ml	1.0%SDS
Igepal	0.25ml	0.5% Igepal
Tween-20	0.25ml	0.5%Tween-20
20mg/ml RNase A	1.0ml	400µg/ml RNase A
(Stock made up in 1mM TrizmaHCl/15mM NaCl – store at -20°C)		
HPLC Water	43.95ml	
Store at 4°C		

Proteinase K

Stock concentration	Final concentration
20mg/ml Proteinase K (>20U/mg) made up in 10mM Tris pH8.0, 40% glycerol – store at -20°C - Invitrogen)	1.2mg/ml

Protocol

1. Add 100µl GenLysis buffer to each well of the plated PBS washed ES cells and then add 3µl Proteinase K to each well
2. Shake for 30sec using a plate shaker.
3. Centrifuge up to 50xg and stop
4. Seal the plate with a sticky seal and incubate at 58°C from 3 hrs to overnight in a sandwich box containing wet towels.
5. Remove the plates from the box and shake for 1min using a plate shaker.
6. Add 5µl RNaseA (20mg/ml), shake for 30sec using a plate shaker and incubate for 30min at 37°C.
7. Add 5µl 5M NaCl.
8. Shake for 30sec using a plate shaker.
9. Transfer 70µl to a 96 well round bottom polypropylene plate
10. Add 80µl 100% Isopropanol.
11. Centrifuge for 1hr at 2,800xg and 4°C and then tip off isopropanol
12. Add 150µl 70% ethanol.
13. Centrifuge 5min at 2,800xg and 4°C
14. Tip off ethanol.
15. Repeat steps 15 -17
16. With plate still inverted, place on a Whatman blotting paper and centrifuge up to 30xg (~15sec) and stop.
17. Add 50µl 0.1xTE and seal.
18. Shake for 30sec, using a plate shaker.
19. Incubate for 30 min at 58°C.

20. Shake for 30sec, using a plate shaker.
21. Leave overnight on the bench
22. For QC use 4 μ l on a 1% agarose gel.
23. Seal with a foil seal and store at -20°C.